


STONE EXTRACTOR

Designed for handling and removal of stone in the urinary tract.
Tipless design allows the capture of the stone at the base of the calyx
reducing the likelihood of tissue damage.





Nitinol and Stainless Steel Basket
2,4 ; 3 and 4 fr.
90 and 120 cm.
With Standard tip and Tipless
Flatwire and Helical
Polymide and Polyethylene tubes



STONE EXTRACTOR



Diameter and length	Basket made of		Format
	Nitinol	Stainless steel	
2,4 fr x 90 cm	NT-24F90R-PI	SS-24F90R-PI	Straight
3,0 fr x 90 cm	NT-3F90R-PI	SS-3F90R-PI	
4,0 fr x 90 cm	NT-4F90R-PI	SS-4F90R-PI	
2,4 fr x 120 cm	NT-24F120R-PI	SS-24F120R-PI	
3,0 fr x 120 cm	NT-3F120R-PI	SS-3F120R-PI	
4,0 fr x 120 cm	NT-4F120R-PI	SS-4F120R-PI	
2,4 fr x 90 cm	NT-24F90H-PI	SS-24F90H-PI	Helical
3,0 fr x 90 cm	NT-3F90H-PI	SS-3F90H-PI	
4,0 fr x 90 cm	NT-4F90H-PI	SS-4F90H-PI	
2,4 fr x 120 cm	NT-24F120H-PI	SS-24F120H-PI	
3,0 fr x 120 cm	NT-3F120H-PI	SS-3F120H-PI	
4,0 fr x 120 cm	NT-4F120H-PI	SS-4F120H-PI	
3,0 fr x 90 cm	NT-3F90R-PE	SS-3F90R-PE	Straight
4,0 fr x 90 cm	NT-4F90R-PE	SS-4F90R-PE	
3,0 fr x 120 cm	NT-3F120R-PE	SS-3F120R-PE	
4,0 fr x 120 cm	NT-4F120R-PE	SS-4F120R-PE	
3,0 fr x 90 cm	NT-3F90H-PE	SS-3F90H-PE	Helical
4,0 fr x 90 cm	NT-4F90H-PE	SS-4F90H-PE	
3,0 fr x 120 cm	NT-3F120H-PE	SS-3F120H-PE	
4,0 fr x 120 cm	NT-4F120H-PE	SS-4F120H-PE	
2,4 fr x 90 cm	SP-24F90R-PI		Straight
2,4 fr x 120 cm	SP-24F120R-PI		
3,0 fr x 90 cm	SP-3F90H-PI		Helical / Tip cero
4,0 fr x 90 cm	SP-4F90H-PI		
3,0 fr x 120 cm	SP-3F120H-PI		
4,0 fr x 120 cm	SP-4F120H-PI		
3,0 fr x 90 cm	SP-3F90H-PE		Helical / Tip cero
4,0 fr x 90 cm	SP-4F90H-PE		
3,0 fr x 120 cm	SP-3F120H-PE		
4,0 fr x 120 cm	SP-4F120H-PE		
3,0 fr x 90 cm	TF-3F90-PI		3-Nail Grasping Forceps
4,0 fr x 90 cm	TF-4F90-PI		
3,0 fr x 120 cm	TF-3F120-PI		
4,0 fr x 120 cm	TF-3F120-PI		
3,0 fr x 90 cm	TF-3F90-PE		3-Nail Grasping Forceps
4,0 fr x 90 cm	TF-4F90-PE		
3,0 fr x 120 cm	TF-3F120-PE		
4,0 fr x 120 cm	TF-3F120-PE		
8,0 fr X 38 cm	GP 8/38		Percutaneous extractor



PI Polyimide sheath
PE Polyethylene sheath