

- ULTRASOFT SURFACE
- WIDE LUMEN
- GREAT FLEXIBILITY
- RESISTEANCE TO ENCRUSTATION
- LOW FRICTION SURFACE

# COMPLIANCE

Double-J ureteral kit



# Double-J ureteral kit


For internal drainage, as from the ureteral pelvic junction to the bladder

- Single use.
- Constructed from a special biocompatible polymer that minimizes sensitivity reactions, allowing an adequate tolerance by patient
- Ultrasoft surface catheter with uniform dimensions.
- Wide lumen.
- Great flexibility designed to avoid collapses.
- Optimum orifice sizing, facilitating a better drainage of fluids.
- Biocompatible markings to indicate depth.

## Components:

- 1 - Double J Catheter.
- 1 - Guidewire.
- 1 - Clamp.
- 1 - Pusher Catheter.

CONFIDANCE



	Open Tip Hydrophilic Catheter	Closed Tip Catheter	Diameter fr.	Length cm / inch.	Guide Diam. inch.
	DAHB - 4520	DCHB- 4520	4,5	20 / 7,87	0,035
	DAHB - 0520	DCHB- 0520	5		0,035
	DAHB - 0620	DCHB- 0620	6		0,038
	DAHB - 0522	DCHB- 0522	5	22 / 8,66	0,035
	DAHB - 0622	DCHB- 0622	6		0,038
	DAHB - 0722	DCHB- 0722	7		0,038
	DAHB - 0822	DCHB- 0822	8		0,038
	DAHB - 0524	DCHB- 0524	5	24 / 9,45	0,035
	DAHB - 0624	DCHB- 0624	6		0,038
	DAHB - 0724	DCHB- 0724	7		0,038
	DAHB - 0824	DCHB- 0824	8		0,038
	DAHB - 0526	DCHB- 0526	5	26 / 10,23	0,035
	DAHB - 0626	DCHB- 0626	6		0,038
	DAHB - 0726	DCHB- 0726	7		0,038
	DAHB - 0826	DCHB- 0826	8		0,038
	DAHB - 0528	DCHB- 0528	5	28 / 11,02	0,035
	DAHB - 0628	DCHB- 0628	6		0,038
	DAHB - 0728	DCHB- 0728	7		0,038
	DAHB - 0828	DCHB- 0828	8		0,038
Multi Length	SCHML-050820		5	8 - 20 cm /	0,035
	SCHML-060820		6	3,15"- 7,87"	0,038
	SCHML-062232		6	20 - 32 cm /	0,038
	SCHML-072232		7	7,87"- 12,6"	0,038

Include the letter "H" to request hydrophilic guidewire.  
Example: DAHBH-0722 ( Catheter and Guidewire Hydrophilic )