

Fadicat

Implantable reservoirs for oncology



- Silicon catheter
- Port having round design and shape to facilitate its placement under the skin with three orifices for suture
- Constructed from polyacetal (Delrin), that guarantees its biocompatibility and enhanced comfort for the patient with the minimum weight
- Designed to allow repeated access.

Use:

- Access in chemotherapy
- Drug administration
- Blood sampling
- Parenteral nutrition
- Transfusions
- Antibiotic therapy
- Long-term treatment of pain
- Chemotherapeutic treatment of hepatic metastases

SPECIFICATIONS

Part N°	Chamber	Silicon catheter	
		Length (cm)	Diameter (fr)
AD09	Adult	75	9
AD09K	Adult	75	9
AD07	Adult	75	7
AD07K	Adult	75	7
AD05	Pediatric	40	5
AD05K	Pediatric	40	5

Part N°: AD09 / AD07 / AD05 Reservoir And Catheter Only.

KIT	AD09K	AD07K	AD05K
Puncture Needle	18 Ga	18 Ga	18 Ga
Guide wire	0,035" x 45 cm	0,035" x 45 cm	0,035" x 45 cm
Pealable Introducer	10 fr.	8 fr.	6 fr.
Tunneler	25 cm	25 cm	25 cm
5 cc Syringe	1	1	1
10 cc Syringe	1	1	1
Scalpel Sheet	1	1	1
Huber Needle	1" x 21 Ga	1" x 21 Ga	1" x 23 Ga
Hypodermic Needle	40/8	40/8	25/6

Fadicat



Manufactured by : KFF s.a.
Espora 49 - Sarandí (C.P. 1872)
Buenos Aires - Argentina
Phone: (54-11) 4265-4515
E-mail: info@kffmed.com
Home page: www.kffmed.com