

Gastrocat

Set for radiological
percutaneous gastrostomy



- Low profile catheter.
- Easier insertion.
- Larger lumen.
- Retention system.
- Extended durability.

Gastrocat

Set for radiological percutaneous gastrostomy

Feasibility of collocation is about 100 % of the cases.

Components:

- Gastrostomy catheter.
- Rigid Cannula.
- Needle to insert anchorage points.
- Set of three percutaneous points.
- Introducer needle, 18 G x 7 cm.
- J-tip guide wire.
- Set of tissue dilators.
- Disk to anchor catheter.
- 60 cc syringe Toomey type.

Code		Catheter diameter	Use
SG - 18	Set	18 Fr.	Adult
CG - 18	Catheter		
SG - 14	Set	14 Fr.	
CG - 14	Catheter		
SG - 12	Set	12 Fr.	
CG - 12	Catheter		
SGP - 14	Set	14 Fr.	Pediatric
CGP - 14	Catheter		
SGP - 12	Set	12 Fr.	
CGP - 12	Catheter		