

# Gastrotop

## Gastrostomy Button

Fully constructed from biocompatible silicone.

Low profile.

Antireflux valve.

Easy placement.

Easy to handle for patients.



# Gastrotop

## Gastrostomy Button

- An effective therapeutic alternative with low placement complexity.
- Designed for long-term enteral feeding for patients who have difficulties in receiving oral nutrition.
- Completely manufactured with 100% medical grade silicone.
- With a one-way valve that prevents reflux and allows the flow of liquids and semisolids.
- Extragastric foldable sector with integrated seal.
- Better comfort for the patient.
- Easy-to-use.
- Fewer hospitalization days.
- Home care with minimum assistance.
- Low-profile device for a better aesthetic appearance and comfort.



| Part N° |       | Diámetro |         |         |
|---------|-------|----------|---------|---------|
|         |       | 18 Fr.   | 20 Fr.  | 24 Fr.  |
| Length  | 12 mm | GT-1812  |         |         |
|         | 17 mm | GT-1817  | GT-2017 | GT-2417 |
|         | 20 mm |          | GT-2020 |         |
|         | 24 mm | GT-1824  | GT-2024 | GT-2424 |
|         | 28 mm | GT-1828  | GT-2028 | GT-2428 |
|         | 34 mm | GT-1834  | GT-2034 | GT-2434 |
|         | 44 mm |          | GT-2044 | GT-2444 |

| Part N° |  | Feeding Tube |           |           |
|---------|--|--------------|-----------|-----------|
|         |  | 18 Fr.       | 20 Fr.    | 24 Fr.    |
|         |  | TAGBGE-18    | TAGBGE-20 | TAGBGE-24 |

| Part N° |  | Descompression Tube |             |             |
|---------|--|---------------------|-------------|-------------|
|         |  | 18 Fr.              | 20 Fr.      | 24 Fr.      |
|         |  | TDBGGE-18xx         | TDBGGE-20xx | TDBGGE-24xx |

xx : button length.  
Ej.: a GT-1824 GASTROTOP button,  
requires a TDBGGE-1824 descompressor.

Components : Feeding tube- Descompression tube (for gastric descompression) - Feeding syringe - Indicator - Placement extension set