

GUARDIAN

Semi-implantable central
venous catheter

2



GUARDIAN

Semi-implantable central venous catheter

2

The silicone central venous catheter is designed for chronic patients in need of repeated access to central venous systems for medicine, nutritional solutions, blood products and blood sample taking.

Available in one, two and three lumens, and different diameter .

Catheter is made of medical grade silicone to offer state-of-the-art bio-compatibility and excellent adaptation to the vessel.

Totally radiopaque allows optimal view during implantation and makes possible to verify correct positioning from the radiological perspective.

A polyester cuff allows tissue growth to reduce bacteria migration; help for further infections and prevents accidental removal of catheter.

Part N°	N° Lumen	Diameter FR	Length mm	Model
HICKES-4290	1	4,2	90	Unilum H
HICKES-0590	1	5	90	
HICKES-0690	1	6	90	
HICKES-0790	1	7	90	
HICKES-9590	1	9,5	90	
HICKED-4290	2	4,2	90	Bilum H
HICKED-0590	2	5	90	
HICKED-0790	2	7	90	
HICKED-9590	2	9,5	90	
HICKED-1190	2	11	90	
HICKED-1490	2	14	90	Trilum H
HICKET-1290	3	12,5	90	

Components.

1 Catheter with polyester cuff
1 Guidewire
1 Tunneler (6 to 14 Fr.)
1 Introducer catheter on needle (4, 2 and 5 Fr.)
1 Peelable introducer with dilator
1 Introducer needle
1 Scalpel
Injection caps
Syringe with hypodermic needle
Clamps

The logo for KFF s.a. features the letters 'KFF' in a bold, stylized font. The 'K' is formed by two vertical bars connected at the top, and the 'FF' consists of two vertical bars with horizontal bars extending from their tops. To the right of the 'KFF' is the text 's.a.' in a smaller, sans-serif font.