


Only for Brazil

KFF s.a.

HIDROGUIDE

Hydrophilic coated



Steerability
Pushability

Torque control 1:1

A GUIDEWIRE WITH EASY NAVEGABILITY

| Part N° | Diameter | Length | Tip | Tip length |
|-------------|----------|--------|----------|------------|
| RHHDE-32150 | 0.032" | 150 cm | Straight | 3 cm |
| AHHDE-32150 | 0.032" | 150 cm | Angle | 3 cm |
| RHHDE-35150 | 0.035" | 150 cm | Straight | 3 cm |
| AHHDE-35150 | 0.035" | 150 cm | Angle | 3 cm |
| RHHDE-35180 | 0.035" | 180 cm | Straight | 3 cm |
| AHHDE-35180 | 0.035" | 180 cm | Angle | 3 cm |