


LONG TERM SILICONE CATHETER

*INNOVA*



An effective solution  
for hemodialysis

**KFF** s.a.

## LONG TERM SILICONE CATHETER

# INNOVA

BIOCOMPATIBLE MATERIAL, OPTIMAL COMFORT FOR THE PATIENT AND GOOD ADAPTATION TO THE VASCULAR ANATOMY.

ROTATING SUTURE WING.

INDICATION OF PRIMING VOLUME.

POLYESTER CUFF FIXATION.

DISTAL HOLES LOCATIONS ACHIEVES RECIRCULATION BLOOD REDUCTION.

Code N°	Diameter fr.	Length cm.	Catheter Type
BLDS-0810	8	10	Straight
BLDS-0818	8	18	Straight
BLDS-0824	8	24	Straight
BLDS-12528	12,5	28	Straight
BLDS-12532	12,5	32	Straight
BLDS-12528-C	12,5	28	Curved
BLDS-12532-C	12,5	32	Curved
BLDS-14524	14,5	24	Straight
BLDS-14528	14,5	28	Straight
BLDS-14532	14,5	32	Straight

### Components:

- 1 Double lumen silicone catheter with cuff.
- 1 Introducer needle.
- 1 Metallic guide wire.
- 2 Vessel dilators.
- 1 Peelable introducer sheath.
- 1 Tunneler with sleeve.
- 2 Injection caps.
- 1 Syringe with hypodermic needle.
- 1 Semipermeable dressing.
- 1 Scalpel.