


INNOVA

For safe, effective treatment

New line of two-lumens
catheters for
hemodialysis and
plasmapheresis



KFF s.a.

New line of two-lumens catheters for hemodialysis and plasmapheresis

INNOVA

For safe, effective treatment

Made of biocompatible polyurethane.

High flow: The INNOVA catheter line allows a high flow, even when working at low pressures, and thus maximizes treatment effectiveness.

For subclavian or jugular vein access.

Straight or curved branches.

The curved branches option, ideal for internal jugular vein access, make the catheter more comfortable for the patient.

Rotating suture wing.

Depth indicator.

Size indication (gauge and length).



Kit Components
1 Double-lumen catheter, 12 fr.
1 Introducer needle, 18 G x 7 cm
1 Guide wire, 0,038" x 70 cm
1 Dilator
2 Injection caps
1 Scalpel
1 Syringe with hypodermic needle

Part N°	Diameter	Useful length	Extensions	Tip	Configuration
KIHR-1215B	12 Fr.	15 cm	Straight	Std	Kit
KIHR-1216B	12 Fr.	16 cm	Straight	Std	Kit
KIHR-1220B	12 Fr.	20 cm	Straight	Std	Kit
KIHC-1215B	12 Fr.	15 cm	Curved	Std	Kit
KIHC-1216B	12 Fr.	16 cm	Curved	Std	Kit
KIHC-1220B	12 Fr.	20 cm	Curved	Std	Kit
KIHR-1215SB	12 Fr.	15 cm	Straight	Soft	Kit
KIHR-1216SB	12 Fr.	16 cm	Straight	Soft	Kit
KIHR-1220SB	12 Fr.	20 cm	Straight	Soft	Kit
KIHC-1215SB	12 Fr.	15 cm	Curved	Soft	Kit
KIHC-1216SB	12 Fr.	16 cm	Curved	Soft	Kit
KIHC-1220SB	12 Fr.	20 cm	Curved	Soft	Kit
IHR-1215	12 Fr.	15 cm	Straight	Std	Catheter
IHR-1216	12 Fr.	16 cm	Straight	Std	Catheter
IHR-1220	12 Fr.	20 cm	Straight	Std	Catheter
IHC-1215	12 Fr.	15 cm	Curved	Std	Catheter
IHC-1216	12 Fr.	16 cm	Curved	Std	Catheter
IHC-1220	12 Fr.	20 cm	Curved	Std	Catheter
IHR-1215S	12 Fr.	15 cm	Straight	Soft	Catheter
IHR-1216S	12 Fr.	16 cm	Straight	Soft	Catheter
IHR-1220S	12 Fr.	20 cm	Straight	Soft	Catheter
IHC-1215S	12 Fr.	15 cm	Curved	Soft	Catheter
IHC-1216S	12 Fr.	16 cm	Curved	Soft	Catheter
IHC-1220S	12 Fr.	20 cm	Curved	Soft	Catheter

FOR VALVED SYRINGE ADD LETTER "V" AT THE END OF THE PART N°

EBRO-INN-001
NOV 2014