

CATHETER SHEATH INTRODUCERS

Constructed from radio-opaque and biocompatible polymers.
High kink resistance.
Excellent transition between sheath and dilator.



INTRODUCER SET:

- Introducer sheath.
- Haemostatic valve.
- Lateral derivation with three-way stopcock.
- J/straight Guide wire.
- Dilator.
- Introducer needle with appropriate length and caliber.
- Protective sleeve-bearing the letter "F" in its code.
- Radial introducer-bearing the letter "R" in its code.

| Catalogue N° | Measure | Guide diameter | Useful length |
|--------------|---------|----------------|---------------|
| ICB-0510 | 5.0 Fr | .035" | 10 cm |
| ICB-0610 | 6.0 Fr | .035" | 10 cm |
| ICB-0710 | 7.0 Fr | .035" | 10 cm |
| ICB-0810 | 8.0 Fr | .035" | 10 cm |
| ICB-8510 | 8.5 Fr | .035" | 10 cm |
| ICB-0910 | 9.0 Fr | .035" | 10 cm |
| ICB-0515 | 5.0 Fr | .035" | 15 cm |
| ICB-0615 | 6.0 Fr | .035" | 15 cm |
| ICB-0715 | 7.0 Fr | .035" | 15 cm |
| ICB-0815 | 8.0 Fr | .035" | 15 cm |
| ICB-8515 | 8.5 Fr | .035" | 15 cm |
| ICB-0915 | 9.0 Fr | .035" | 15 cm |
| ICB-0523 | 5.0 Fr | .035" | 23 cm |
| ICB-0623 | 6.0 Fr | .035" | 23 cm |
| ICB-0523 R | 5.0 Fr | .021" | 23 cm |
| ICB-0623 R | 6.0 Fr | .021" | 23 cm |
| ICB-0723 | 7.0 Fr | .035" | 23 cm |
| ICB-0723 R | 7.0 Fr | .021" | 23 cm |
| ICB-0823 | 8.0 Fr | .035" | 23 cm |
| ICB-0510 F | 5.0 Fr | .035" | 10 cm |
| ICB-0610 F | 6.0 Fr | .035" | 10 cm |
| ICB-0710 F | 7.0 Fr | .035" | 10 cm |
| ICB-0810 F | 8.0 Fr | .035" | 10 cm |
| ICB-8510 F | 8.5 Fr | .035" | 10 cm |
| ICB-0910 F | 9.0 Fr | .035" | 10 cm |

CATHETER SHEATH INTRODUCERS

Constructed from radio-opaque and biocompatible polymers. High kink resistance. Excellent transition between sheath and dilator.

PEELABLE INTRODUCER KIT:

- Introducer
- J/straight Guide wire
- Dilator.
- Introducer needle

| Catalogue N° | Measure | Guide diameter | Useful length |
|--------------|---------|----------------|---------------|
| IPB - 0810 | 8.0 Fr | .035" | 10 cm |
| IPB - 0910 | 9.0 Fr | .035" | 10 cm |
| IPB - 1010 | 10.0 Fr | .035" | 10 cm |
| IPB - 1110 | 11.0 Fr | .035" | 10 cm |
| IPB - 1210 | 12.0 Fr | .035" | 10 cm |
| IPB - 0815 | 8.0 Fr | .035" | 15 cm |
| IPB - 0915 | 9.0 Fr | .035" | 15 cm |
| IPB - 1015 | 10.0 Fr | .035" | 15 cm |
| IPB - 1115 | 11.0 Fr | .035" | 15 cm |
| IPB - 1215 | 12.0 Fr | .035" | 15 cm |