

# GUARDIAN

## PICC Catheter



**KFF** s.a.

# Guardian

## PICC Catheter

**Designed for antibiotic therapy, drugs for treatment of pain, chemotherapy and parenteral nutrition.**

- Ultra-soft silicone surface guarantees reliability, less discomfort to the patient and better tissue care.
- Excellent bio&hemocompatibility allowing greater permanence.
- Easy insertion.
- Suture wing .
- Depths marks.

Being aware of different medical requirements, 1 or 2 lumens catheters are supplied in different sizes, as well as various insertion systems.

| Part. N°   | Diameter | Usef. Leng. | N° Lumen |
|------------|----------|-------------|----------|
| CCIPU-0230 | 2 Fr.    | 30          | 1        |
| CCIPU-0360 | 3 Fr.    | 60          | 1        |
| CCIPD-0360 | 3 Fr.    | 60          | 2        |
| CCIPU-0460 | 4 fr.    | 60          | 1        |
| CCIPD-0460 | 4 Fr.    | 60          | 2        |
| CCIPU-0560 | 5 Fr.    | 60          | 1        |
| CCIPD-0560 | 5 Fr.    | 60          | 2        |
| CCIPU-0660 | 6 Fr.    | 60          | 1        |
| CCIPD-0660 | 6 Fr.    | 60          | 2        |
| UCCIP-0360 | 3 Fr.    | 60          | 1        |
| DCCIP-0360 | 3 Fr.    | 60          | 2        |
| UCCIP-0460 | 4 Fr.    | 60          | 1        |
| DCCIP-0460 | 4 Fr.    | 60          | 2        |
| UCCIP-0560 | 5 Fr.    | 60          | 1        |
| DCCIP-0560 | 5 Fr.    | 60          | 2        |
| UCCIP-0660 | 6 Fr.    | 60          | 1        |
| DCCIP-0660 | 6 Fr.    | 60          | 2        |



### CCIPU Components

- 1 Single lumen PICC catheter.
- 1 Mandril guide.
- 1 Peelable needle.
- 1 Lateral derivation for irrigation.
- 1 Injection cap.
- 1 Scalpel sheet.
- 1 Syringe with hypodermic needle.
- 1 Semipermeable dressing.

### CCIPD Components

- 1 Double lumen PICC catheter.
- 1 Mandril guide.
- 1 Peelable needle.
- 1 Lateral derivation for irrigation.
- 2 Injection cap.
- 1 Scalpel sheet.
- 1 Syringe with hypodermic needle.
- 1 Semipermeable dressing.

### UCCIP Component (Micro Punction)

- 1 Single lumen PICC catheter.
- 1 45 cm. guidewire (straight / J).
- 1 Introducer needle.
- 1 Peelable introducer with dilator.
- 1 Injection cap.
- 1 Scalpel.
- 1 Syringe with hypodermic needle.
- 1 Semipermeable dressing.

### DCCIP Component (Micro Punction)

- 1 Double lumen PICC catheter.
- 1 45 cm. guidewire (straight / J).
- 1 Injection caps.
- 1 Peelable introducer with dilator.
- 2 Injection caps.
- 1 Scalpel.
- 1 Syringe with hypodermic needle.
- 1 Semipermeable dressing.