



## CENTRAL VENOUS CATHETERIZATION SET

CONSTRUCTED FROM RADIO-OPAQUE  
AND BIOCOMPATIBLE POLYURETHANE.  
SOFT TIP. LOW FRICTION SURFACE.  
DEPTH MARKINGS.

**New!**

0,032" Guidewire  
(Optional)

New Bevel Needle  
More sliding surface  
Better insertion



## CENTRAL VENOUS CATHETERIZATION SET

CONSTRUCTED FROM RADIO-OPAQUE AND BIOCOMPATIBLE POLYURETHANE. SOFT TIP. LOW FRICTION SURFACE. DEPTH MARKINGS.

Model	Part N°	Measure	Useful length	N° of lumen	Contents
<b>UNILUM S</b>	USBA-14G20	14 G	20 cm	1	Catheter with integral suture wing. J/straight Guide wire Vessel dilator/s Introducer needle Syringe 5 cc and hypodermic needle Movable Wing
	USBA-14G30	14 G	30 cm	1	
	USBA-16G20	16 G	20 cm	1	
	USBA-16G30	16 G	30 cm	1	
	USBA-18G13	18 G	13 cm	1	
	USBA-18G20	18 G	20 cm	1	
	USBA-18G30	18 G	30 cm	1	
<b>BILUM S</b>	DSB-0405	4.0 F	5 cm	2	Catheter with integral suture wing. J/straight Guide wire Vessel dilators Introducer needle Syringe 5 cc and hypodermic needle Movable Wing Clamps
	DSB-0408	4.0 F	8 cm	2	
	DSB-0413	4.0 F	13 cm	2	
	DSB-0430	4.0 F	30 cm	2	
	DSB-0505	5.0 F	5 cm	2	
	DSB-0508	5.0 F	8 cm	2	
	DSB-0513	5.0 F	13 cm	2	
	DSB-0520	5.0 F	20 cm	2	
	DSB-0530	5.0 F	30 cm	2	
	DSB-0720	7.0 F	20 cm	2	
	DSB-0730	7.0 F	30 cm	2	
<b>TRILUM S</b>	TSB-5513	5.5 F	13 cm	3	
	TSB-5530	5.5 F	30 cm	3	
	TSB-0720	7.0 F	20 cm	3	
	TSB-0730	7.0 F	30 cm	3	