

CATHETER SHEATH INTRODUCERS

Constructed from radio-opaque and biocompatible polymers.
High kink resistance.
Excellent transition between sheath and dilator.



INTRODUCER SET:

- Introducer sheath.
- Haemostatic valve.
- Lateral derivation with three-way stopcock.
- J/straight Guide wire.
- Dilator.
- Introducer needle with appropriate length and caliber.
- Protective sleeve-bearing the letter "F" in its code.
- Radial introducer-bearing the letter "R" in its code.

Catalogue N°	Measure	Guide diameter	Useful length
ICB-0510	5.0 Fr	.035"	10 cm
ICB-0610	6.0 Fr	.035"	10 cm
ICB-0710	7.0 Fr	.035"	10 cm
ICB-0810	8.0 Fr	.035"	10 cm
ICB-8510	8.5 Fr	.035"	10 cm
ICB-0910	9.0 Fr	.035"	10 cm
ICB-0515	5.0 Fr	.035"	15 cm
ICB-0615	6.0 Fr	.035"	15 cm
ICB-0715	7.0 Fr	.035"	15 cm
ICB-0815	8.0 Fr	.035"	15 cm
ICB-8515	8.5 Fr	.035"	15 cm
ICB-0915	9.0 Fr	.035"	15 cm
ICB-0523	5.0 Fr	.035"	23 cm
ICB-0623	6.0 Fr	.035"	23 cm
ICB-0523 R	5.0 Fr	.021"	23 cm
ICB-0623 R	6.0 Fr	.021"	23 cm
ICB-0723	7.0 Fr	.035"	23 cm
ICB-0723 R	7.0 Fr	.021"	23 cm
ICB-0823	8.0 Fr	.035"	23 cm
ICB-0510 F	5.0 Fr	.035"	10 cm
ICB-0610 F	6.0 Fr	.035"	10 cm
ICB-0710 F	7.0 Fr	.035"	10 cm
ICB-0810 F	8.0 Fr	.035"	10 cm
ICB-8510 F	8.5 Fr	.035"	10 cm
ICB-0910 F	9.0 Fr	.035"	10 cm

CATHETER SHEATH INTRODUCERS

Constructed from radio-opaque and biocompatible polymers. High kink resistance. Excellent transition between sheath and dilator.

PEELABLE INTRODUCER KIT:

- Introducer
- J/straight Guide wire
- Dilator.
- Introducer needle

Catalogue N°	Measure	Guide diameter	Useful length
IPB - 0810	8.0 Fr	.035"	10 cm
IPB - 0910	9.0 Fr	.035"	10 cm
IPB - 1010	10.0 Fr	.035"	10 cm
IPB - 1110	11.0 Fr	.035"	10 cm
IPB - 1210	12.0 Fr	.035"	10 cm
IPB - 0815	8.0 Fr	.035"	15 cm
IPB - 0915	9.0 Fr	.035"	15 cm
IPB - 1015	10.0 Fr	.035"	15 cm
IPB - 1115	11.0 Fr	.035"	15 cm
IPB - 1215	12.0 Fr	.035"	15 cm