

Gastrocat

Set for Radiological
Percutaneous Gastrostomy

- Low Diameter Pigtail Catheter.
- Large Lumen.
- Non-Inflatable Retention System.
- Easier Insertion.
- Greater Durability.



Gastroplexy Points

Gastrocat

Set for Radiological Percutaneous Gastrostomy



CODE	DIAMETER		
CG-12	12 fr.	Catheter	Pigtail Central
CG-14	14 fr.		
CG-18	18 fr.		
SG-12	12 fr.	Set	Pigtail Central
SG-14	14 fr.		
SG-18	18 fr.		
CGP-12	12 fr.	Catheter	Pigtail Distal
CGP-14	14 fr.		
CGP-18	18 fr.		
SGP-12	12 fr.	Set	Pigtail Distal
SGP-14	14 fr.		
SGP-18	18 fr.		

Componentes:

- Gastrostomy catheter.
- Rigid Cannula.
- Needle to insert anchorage points.
- Set of three percutaneous points.
- Introducer needle, 18 G x 7 cm.
- J-tip guidewire.
- Set of tissue dilators.
- Disk to anchor catheter.
- 60 cc syringe Toomey type.

